

## Investigation 6: Spectrophotometric Measurement of an Equilibrium Constant

**Focus Question:** How can we determine the equilibrium constant for bromothymol blue using spectrophotometric data?

### Pre-lab reading

[Chemistry: an Atoms-Focused Approach:](#) Sections 14.1-14.2, 14.7

Primers:

[Spectrophotometry](#)

[SpectraVis Plus spectrophotometer](#)

[Volumetric glassware use – General](#)

[Volumetric glassware use – volumetric flask](#)


[Volumetric glassware use – volumetric pipet](#)

### Safety and Waste Disposal

- HCl and NaOH solutions are a strong acid and base, caution should be used when handling.
- All solutions used in this lab can be disposed of down the drain with lots of water.

### Background

Bromothymol blue is an acid-base indicator that is yellow at pH values less than 6 and blue at pH values above 8. Between pH values of 6 and 8, bromothymol blue is some shade of green. It is more yellow-green near pH 6 and more blue-green near pH 8. Its chemical structure is shown in Figure 1.




**Figure 1:** Chemical structure of bromothymol blue (3', 3''-dibromothymolsulfonephthalein) as it occurs under acidic conditions. Its chemical formula is  $\text{HC}_{27}\text{H}_{27}\text{Br}_2\text{O}_5\text{S}$ .

Like all acid-base indicators, bromothymol blue is a weak acid as shown in the following equilibrium in which HBB represents the protonated (yellow) form of bromothymol blue and  $\text{BB}^-$  represents the blue conjugate base form.



Absorbance spectra (plots of absorbance versus wavelength) for the yellow and blue forms are shown in Figure 2. The blue form,  $\text{BB}^-$ , has an absorbance maximum at about 616 nm. The yellow form, HBB, has its maximum absorbance at 432 nm. In this experiment, we will measure the absorbance of the yellow form at 453 nm, where the absorbance is still strong and the absorbance of the blue  $\text{BB}^-$  is minimal (as can be seen in the spectra in Figure 2 at 432 nm,  $\text{BB}^-$  also absorbs).



**Figure 2:** Visible absorbance spectra of bromothymol blue at low pH (yellow, solid line) and high pH (blue, dashed line).

The equilibrium expression,  $K_c$ , for this equilibrium (Eq. 1) is:

$$K_c = \frac{[BB^-][H_3O^+]}{[HBB]} \quad (2)$$

The value of  $K_c$  should be independent of all factors except a change in temperature. At high pH, the concentration of the blue form,  $[BB^-]$ , is large and  $[HBB]$  is small, and at low pH,  $[HBB]$  is large and  $[BB^-]$  is small. In this experiment, you will be working with a solution that is green in color, so that neither  $[HBB]$  nor  $[BB^-]$  is much larger than the other.

Both the yellow and blue forms of bromothymol blue have Beer's law expressions ( $A = \epsilon \cdot b \cdot c$ ) that relate absorbance ( $A$ ) to concentration ( $c$  or  $[ ]$ , in molarity) at their respective wavelengths:

$$\text{For the yellow form, HBB: } A_{453nm} = \epsilon_{453nm}^{yellow} \cdot b \cdot [HBB] \quad (3)$$

$$\text{For the blue form, BB: } A_{616nm} = \epsilon_{616nm}^{blue} \cdot b \cdot [BB^-] \quad (4)$$

where  $\epsilon$  is the molar absorptivity ( $L \cdot mol^{-1} \cdot cm^{-1}$ ) and  $b$  is the pathlength in cm. By knowing the molar absorptivity of each form of bromothymol blue and the path length of our cuvette it is possible to use the measured absorbance to calculate the equilibrium concentrations of HBB and  $BB^-$ . However, what is typically done is to make a calibration curve or Beer's law plot at the wavelength of interest using solutions of known concentration. This plot can then be used to determine the concentration of a species in an equilibrium solution.

There are five values that you must obtain to determine  $K_c$  for bromothymol blue equilibrium reaction. These are summarized in Table 1. You will measure three of these values for a solution that is green in color, meaning that both yellow HBB and blue  $BB^-$  are present in reasonable quantities (not almost zero as would be the case in either

the blue or yellow solutions at high and low pH values, respectively).

Table 1: Measured quantities

Quantity	How do we measure it?
$[\text{H}_3\text{O}^+]_{\text{eq}}$	Use a pH meter to determine the pH of a green solution and convert value to concentration.
$A_{453}$	Measure the absorbance of a green solution at 453nm
$A_{616}$	Measure the absorbance of a green solution at 616nm
[HBB]	Beer's Law plot of yellow solution (See below)
[BB <sup>-</sup> ]	Beer's Law plot of a blue solution (See below)

The concentration of the yellow form [HBB], will be determined using a Beer's law plot that is obtained for solutions that are completely yellow (low pH; ~100% HBB and negligible BB<sup>-</sup>) and the concentration of the blue form [BB<sup>-</sup>] will be determined from another Beer's law plot that is obtained from solutions that are completely blue (high pH; ~100% BB<sup>-</sup>).

## Procedure:

### Part 1: Making Beer's Law Plots

1. Obtain enough stock solution ( $4 \times 10^{-4}$  M) of bromothymol blue to make two sets of solutions with concentrations between  $4 \times 10^{-5}$  M and  $4 \times 10^{-6}$  M bromothymol blue. One set (low pH, yellow) will be made using 0.1 M HCl as the diluent and the other set (high pH, blue) will be made using 0.1 M NaOH as the diluent. Use 10 mL volumetric flasks and the appropriate pipets to make four solutions in the correct concentration range using 0.1 M HCl to dilute to the mark. Make a second set of four solutions in 10 mL volumetric flasks in the correct concentration range using 0.1 M NaOH to dilute to the mark.
2. Measure the absorbance of each solution (yellow at 453 nm and blue at 616 nm) and create two Beer's Law plots (absorbance vs. concentration). Fit the data with a linear least squares fit, if any data does not fit the trend re-measure absorbance and remake solutions as needed. The slopes of the best-fit linear equations will be used to obtain a value for the concentrations of HBB and BB<sup>-</sup> in the equilibrium solutions below.

### Part 2: Equilibrium measurements

1. Add 40 drops of bromothymol blue indicator (stock solution) to 10.0 mL of the pH = 7 buffer solution in a small beaker and stir to mix the indicator uniformly through the solution. The solution should appear green in color.
2. If the solution is green, you can skip this step. If the solution is blue, add 1 M HCl(aq), a drop at a time, with thorough mixing, until you obtain a shade of green. If the solution is yellow, add 1 M NaOH (aq), a drop at a time, again with thorough mixing, until you obtain a shade of green.
3. Determine the pH of the green buffer solution using the pH probe. Record the pH value in your notebook.
4. Take a portion out of the beaker and place in a cuvette to measure the absorbance at 453 and 616 nm of this green solution and record the values in your notebook. Return the portion to the sample beaker.
5. Add a few drops of either 1 M NaOH or 1 M HCl to your mixture and repeat steps 3 and 4 until you have six trials for green solutions with different pH's (possibly 3 slightly more basic and 3 slightly more acidic).
6. Use the Beer's law plots to determine the concentration of HBB and BB<sup>-</sup> for each trial. Use the relation  $[\text{H}_3\text{O}^+] = 10^{-\text{pH}}$  to determine the concentration of hydronium ion.
7. Using the values determined calculate the equilibrium constant for the bromothymol blue reaction for each of your trials at different pH. Average the values to obtain a value of  $K_c$  for bromothymol blue reaction.

## References

Atkins, P.; Jones, L. "Chemical Principles: The Quest for Insight", 6th ed.; Freeman: New York. **2013**.